

Comprehensive Hearing Loss and Deafness Panel

Test code: EA0501

Is a 239 gene panel that includes assessment of non-coding variants.

In addition, it also includes the maternally inherited mitochondrial genome.

Is ideal for patients with a clinical suspicion of syndromic or non-syndromic genetic hearing loss.

Hearing loss is a genetically and clinically heterogenous group of diseases and syndromes and may be classified in several different ways. This Panel includes comprehensively genes associated with both syndromic and non-syndromic hearing loss. In addition to protein coding regions, two disease causing intronic variants of *HGF* gene are targeted in this Panel. Inheritance of these disorders may be autosomal recessive or dominant as well as X-linked. This comprehensive Panel includes Waardenburg Syndrome Panel, Pendred Syndrome Panel, Usher Syndrome Panel, Stickler Syndrome Panel, Alport Syndrome Panel, Branchio-Oto-Renal Panel, Syndromic Hearing Loss Panel and Non-Syndromic Hearing Loss Panel.

About Comprehensive Hearing Loss and Deafness

Hearing loss is a genetically very heterogenous group of phenotypes varying in severity and causes. Non-syndromic sensorineural hearing loss is a partial or total loss of hearing that occurs without other associated clinical findings. In syndromic hearing loss, symptoms affecting other parts of the body occur in addition to hearing impairment or deafness. Sensorineural hearing loss can be unilateral or bilateral and it can be stable or progressive. In addition, the loss may appear with various intensity to high, middle or low tones. It is estimated that approximately 60-80% of congenital hearing loss is genetic in origin. Some 60%-to-70% of congenital hereditary hearing impairment have a non-syndromic origin. The prevalence of non-syndromic hearing loss is 3-4:10,000 neonates and increases with age. In many populations, mutations in *GJB2* are the most prevalent explaining up to 50% of all non-syndromic hearing losses. Altogether syndromic hearing loss accounts for 20% to 30% of congenital hearing loss and deafness and the combined prevalence of syndromic hearing loss is approximately 1-2:10,000.

Availability

Results in 3-4 weeks

Gene set description

Genes in the Comprehensive Hearing Loss and Deafness Panel and their clinical significance

Gene	Associated phenotypes	Inheritance	ClinVar	HGMD
ABHD12	Polyneuropathy, hearing loss, ataxia, retinitis pigmentosa, and cataract	AR	16	20
ACTG1*	Deafness, Baraitser-Winter syndrome	AD	27	47
ADCY1	Deafness	AR	1	1
ADGRV1	Febrile seizures, familial, 4, Usher syndrome, type IIC	AR	71	236
AIFM1	Deafness, Combined oxidative phosphorylation deficiency 6, Cowchock syndrome	XL	27	31
ALMS1*	Alström syndrome	AR	197	302
ANKH	Calcium pyrophosphate deposition disease (familial chondrocalcinosis type 2), Craniometaphyseal dysplasia autosomal dominant type	AD	13	20

Blueprint Genetics

ARSG	Usher syndrome, type IV	AR	1	1
ATP2B2	Sensorineural hearing loss	AD	3	7
ATP6V1B1	Renal tubular acidosis with deafness	AR	15	56
ATP6V1B2	Deafness, congenital, with onychodystrophy, autosomal dominant, Zimmermann-Laband syndrome 2	AD	6	3
BCS1L	Bjornstad syndrome, GRACILE syndrome, Leigh syndrome, Mitochondrial complex III deficiency, nuclear type 1	AR	42	37
<u>BDP1*</u>	Hearing loss	AD/AR	1	1
BSND	Sensorineural deafness with mild renal dysfunction, Bartter syndrome	AR	10	20
BTD	Biotinidase deficiency	AR	170	247
C10ORF2	Perrault syndrome, Mitochondrial DNA depletion syndrome, Progressive external ophthalmoplegia with mitochondrial DNA deletions, autosomal dominant, 3	AD/AR	37	80
CABP2	Deafness	AR	1	6
CACNA1D	Primary aldosteronism, seizures, and neurologic abnormalities, Sinoatrial node dysfunction and deafness	AD/AR	7	8
CATSPER2	Male Infertility, Deafness	AR	2	7
CCDC50	Deafness	AD	1	4
CD151	Raph blood group, Nephropathy with pretibial epidermolysis bullosa and deafness	AR	1	3
CD164	Deafness, autosomal dominant 66	AD	1	1
CDC14A	Deafness, autosomal recessive 105	AR	7	9
CDC42	Takenouchi-Kosaki syndrome, Noonan-syndrome like phenotype	AD	11	9
CDH23	Deafness, Usher syndrome, type 1D	AR	94	358
CDK9		AR		1
CDKN1C	Beckwith-Wiedemann syndrome, IMAGE syndrome	AD	35	81
CEACAM16	Deafness	AD/AR	4	4
CEP250	Cone rod dystrophy and hearing loss	AR		5
CEP78	Cone rod dystrophy and hearing loss	AR	7	9
CHD7	Isolated gonadotropin-releasing hormone deficiency, CHARGE syndrome	AD	276	860
CHSY1	Temtamy preaxial brachydactyly syndrome	AR	6	16
CIB2	Deafness, Usher syndrome type IJ	AR	5	18
CLDN14	Deafness	AR	11	12

Blueprint Genetics

CLIC5	Deafness	AR	1	2
CLPP	Deafness	AR	4	13
CLRN1	Retinitis pigmentosa, Usher syndrome, type 3A	AR	24	39
COCH	Deafness	AD	14	29
COL11A1	Marshall syndrome, Fibrochondrogenesis, Stickler syndrome type 2	AD/AR	34	94
COL11A2	Weissenbacher-Zweymuller syndrome, Deafness, Otospondylomegaepiphyseal dysplasia, Fibrochondrogenesis, Stickler syndrome type 3 (non-ocular)	AD/AR	29	57
COL2A1	Avascular necrosis of femoral head, Rhegmatogenous retinal detachment, Epiphyseal dysplasia, with myopia and deafness, Czech dysplasia, Achondrogenesis type 2, Platyspondylic dysplasia Torrance type, Hypochondrogenesis, Spondyloepiphyseal dysplasia congenital (SEDC), Spondyloepimetaphyseal dysplasia (SEMD) Strudwick type, Kniest dysplasia, Spondyloperipheral dysplasia, Mild SED with premature onset arthrosis, SED with metatarsal shortening, Stickler syndrome type 1	AD	180	561
COL4A3	Alport syndrome, Hematuria, benign familial	AD/AR	123	264
COL4A4	Alport syndrome	AD/AR	110	232
COL4A5	Alport syndrome	XL	704	992
COL4A6	Deafness, with cochlear malformation	XL	11	5
COL9A1	Stickler syndrome recessive type, Multiple epiphyseal dysplasia type 6 (EDM6)	AD/AR	9	6
COL9A2	Stickler syndrome, Multiple epiphyseal dysplasia type 2 (EDM2)	AD/AR	7	12
COL9A3	Multiple epiphyseal dysplasia type 3 (EDM3), Stickler syndrome recessive type	AD/AR	10	14
CRYM	Deafness	AD	2	4
DCAF17	Woodhouse-Sakati syndrome	AR	14	14
DCDC2	Deafness, Nephronophthisis, Sclerosing cholangitis, neonatal	AR	13	9
DFNA5	Deafness	AD	7	13
DFNB31	Deafness, Usher syndrome, type 2D	AR	12	31
DFNB59	Deafness	AR	12	20
DIABLO	Deafness	AD	1	2
DIAPH1	Deafness, Seizures, cortical blindness, and microcephaly syndrome (SCBMS)	AD/AR	10	15
DIAPH3	Non-syndromic sensorineural deafness	AD	1	9
DLX5	Split-hand/foot malformation with sensorineural hearing loss	AR	3	9

Blueprint Genetics

DMXL2	Deafness, autosomal dominant, 71, Polyendocrine-polyneuropathy syndrome, Epileptic encephalopathy, early infantile	AD/AR	2	6
DNMT1	Neuropathy, hereditary sensory, Cerebellar ataxia, deafness, and narcolepsy	AD	9	20
DSPP	Dentin dysplasia, Dentinogenesis imperfecta, Deafness, with dentinogenesis imperfecta	AD	11	53
EDN3	Hirschsprung disease, Central hypoventilation syndrome, congenital, Waardenburg syndrome	AD/AR	7	21
EDNRB	Hirschsprung disease, ABCD syndrome, Waardenburg syndrome	AD/AR	12	66
EIF3F	Intellectual disability	AR		
ELMOD3	Deafness	AR	1	2
EPS8	Deafness	AR	2	2
EPS8L2	Deafness, autosomal recessive 106	AR	2	2
ESPN*	Deafness	AD/AR	12	15
ESRRB	Deafness	AR	12	19
EYA1	Otofaciocervical syndrome, Branchiootic syndrome, Branchiootorenal syndrome	AD	56	218
EYA4	Dilated cardiomyopathy (DCM), Deafness, autosomal dominant 10	AD	15	28
FAM136A	Sensorineural hearing loss	AD	1	2
FAM65B	Deafness	AR	1	2
FDXR	Auditory neuropathy and optic atrophy	AR	5	19
FGF3	Deafness, congenital with inner ear agenesis, microtia, and microdontia	AR	13	20
FGFR3	Lacrimoauriculodentodigital syndrome, Muenke syndrome, Crouzon syndrome with acanthosis nigricans, Camptodactyly, tall stature, and hearing loss (CATSHL) syndrome, Achondroplasia, Hypochondroplasia, Thanatophoric dysplasia type 1, Thanatophoric dysplasia type 2, SADDAN	AD/AR	54	77
FITM2	Dystonia, Deafness	AR		1
FOXI1	Pendred syndrome, Enlarged vestibular aqueduct	AR	1	11
GATA3	Hypomagnesemia, renal	AD	22	86
GIPC3	Deafness	AR	9	20
GJA1*	Oculodentodigital dysplasia mild type, Oculodentodigital dysplasia severe type, Syndactyly type 3	AD/AR	31	107
GJB2	Deafness, Bart-Pumphrey syndrome, Keratoderma, palmoplantar, with deafness, Vohwinkel syndrome, Hystrix-like ichthyosis with deafness, Keratitis-ichthyosis-deafness syndrome	AD/AR/Digenic	133	405

Blueprint Genetics

GJB3	Deafness, Erythrokeratoderma variabilis et progressiva 1, Deafness, autosomal dominant 2B	AD/AR	11	40
GJB6	Deafness, Deafness, autosomal dominant 3B, Ectodermal dysplasia, hidrotic (Clouston syndrome)	AD/AR	10	33
GPSM2	Deafness, Chudley-McCullough syndrome	AR	18	11
GRHL2	Ectodermal dysplasia/short stature syndrome, Deafness, autosomal dominant 28, Corneal dystrophy, posterior polymorphous	AD/AR	12	12
GRXCR1	Deafness	AR	8	9
GRXCR2	Deafness	AR	1	2
HARS*	Charcot-Marie-Tooth disease, axonal, type 2W, Usher syndrome, type 3B	AR	6	12
HARS2	Perrault syndrome	AR	7	3
HGF	Deafness	AR	4	10
HOMER2	Deafness	AD	2	1
HOXB1	Facial paresis, hereditary congenital	AR	3	6
HSD17B4	Perrault syndrome, D-bifunctional protein deficiency	AR	60	99
ILDR1	Deafness	AR	8	27
KARS	Charcot-Marie-Tooth disease	AR	9	23
KCNE1	Long QT syndrome, Jervell and Lange-Nielsen syndrome	AD/AR/Digenic	11	46
KCNJ10	Seizures, sensorineural deafness, ataxia, mental retardation, and electrolyte imbalance (SESAME syndrome), Pendred syndrome, Enlarged vestibular aqueduct	AR/Digenic	13	29
KCNQ1	Short QT syndrome, Long QT syndrome, Atrial fibrillation, Jervell and Lange-Nielsen syndrome	AD/AR/Digenic	298	631
KCNQ4	Deafness	AD	28	37
KIT	Gastrointestinal stromal tumor, Piebaldism	AD	79	116
LARS2	Perrault syndrome, Hydrops, lactic acidosis, and sideroblastic anemia (HLASA)	AR	14	14
LHFPL5	Deafness	AR	7	10
LMX1A	Hearing loss	AD/AR	1	4
LOXHD1	Deafness	AR	26	60
LRP2	Donnai-Barrow syndrome, Faciooculoacousticorenal syndrome	AR	24	38
LRTOMT	Deafness	AR	7	17
MAN2B1	Mannosidosis, alpha B, lysosomal	AR	63	149
MANBA	Mannosidosis, lysosomal	AR	16	19

Blueprint Genetics

MARVELD2	Deafness	AR	9	17
MET	Deafness, Renal cell carcinoma, papillary, Osteofibrous dysplasia, susceptibility to	AD/AR	20	34
MGP	Keutel syndrome	AR	5	8
MIR96	Deafness	AD	2	4
MITF	Tietz albinism-deafness syndrome, Waardenburg syndrome, Coloboma, osteopetrosis, microphthalmia, macrocephaly, albinism, and deafness (COMMAD)	AD/AR	32	58
MPZL2	Sensorineural hearing loss	AR		4
MSRB3	Deafness	AR	5	2
MT-ATP6	Neuropathy, ataxia, and retinitis pigmentosa, Leber hereditary optic neuropathy, Ataxia and polyneuropathy, adult-onset, Cardiomyopathy, infantile hypertrophic, Leigh syndrome, Striatonigral degeneration, infantile, mitochondrial	Mitochondrial	19	
MT-ATP8	Cardiomyopathy, apical hypertrophic, and neuropathy, Cardiomyopathy, infantile hypertrophic	Mitochondrial	4	
MT-CO1	Myoglobinuria, recurrent, Leber hereditary optic neuropathy, Sideroblastic anemia, Cytochrome C oxidase deficiency	Mitochondrial	17	
MT-CO2	Cytochrome c oxidase deficiency	Mitochondrial	8	
MT-CO3	Cytochrome c oxidase deficiency, Leber hereditary optic neuropathy	Mitochondrial	9	
MT-CYB		Mitochondrial	69	
MT-ND1	Mitochondrial myopathy, encephalopathy, lactic acidosis, and stroke-like episodes, Leber hereditary optic neuropathy, Leber optic atrophy and dystonia	Mitochondrial	21	
MT-ND2	Leber hereditary optic neuropathy, Mitochondrial complex I deficiency	Mitochondrial	6	
MT-ND3	Leber optic atrophy and dystonia, Mitochondrial complex I deficiency	Mitochondrial	7	
MT-ND4	Leber hereditary optic neuropathy, Leber optic atrophy and dystonia, Mitochondrial complex I deficiency	Mitochondrial	11	
MT-ND4L	Leber hereditary optic neuropathy	Mitochondrial	2	
MT-ND5	Myoclonic epilepsy with ragged red fibers, Mitochondrial myopathy, encephalopathy, lactic acidosis, and stroke-like episodes, Leber hereditary optic neuropathy, Mitochondrial complex I deficiency	Mitochondrial	19	
MT-ND6	Mitochondrial myopathy, encephalopathy, lactic acidosis, and stroke-like episodes, Oncocytoma, Leber hereditary optic neuropathy, Leber optic atrophy and dystonia, Mitochondrial complex I deficiency	Mitochondrial	16	
MT-RNR1	Deafness, mitochondrial	Mitochondrial	3	
MT-RNR2	Chloramphenicol toxicity/resistance	Mitochondrial	2	
MT-TA		Mitochondrial	4	

Blueprint Genetics

MT-TC	Mitochondrial myopathy, encephalopathy, lactic acidosis, and stroke-like episodes	Mitochondrial	3	
MT-TD		Mitochondrial	1	
MT-TE	Diabetes-deafness syndrome, Mitochondrial myopathy, infantile, transient, Mitochondrial myopathy with diabetes	Mitochondrial	5	
MT-TF	Myoclonic epilepsy with ragged red fibers, Nephropathy, tubulointerstitial, Encephalopathy, mitochondrial, Epilepsy, mitochondrial, Myopathy, mitochondrial, Mitochondrial encephalomyopathy with lactic acidosis and stroke-like episodes	Mitochondrial	7	
MT-TG		Mitochondrial	3	
MT-TH		Mitochondrial	4	
MT-TI		Mitochondrial	7	
MT-TK		Mitochondrial	5	
MT-TL1	Cytochrome c oxidase deficiency, Myoclonic epilepsy with ragged red fibers, Mitochondrial myopathy, encephalopathy, lactic acidosis, and stroke-like episodes, Diabetes-deafness syndrome, Cyclic vomiting syndrome, SIDS, susceptibility to	Mitochondrial	14	
MT-TL2	Mitochondrial multisystemic disorder, Progressive external ophthalmoplegia	Mitochondrial	5	
MT-TM	Leigh syndrome, Mitochondrial multisystemic disorder	Mitochondrial	1	
MT-TN	Progressive external ophthalmoplegia, Mitochondrial multisystemic disorder	Mitochondrial	3	
MT-TP		Mitochondrial	2	
MT-TQ	Mitochondrial multisystemic disorder	Mitochondrial	2	
MT-TR	Encephalopathy, mitochondrial	Mitochondrial	2	
MT-TS1	Myoclonic epilepsy with ragged red fibers, Mitochondrial myopathy, encephalopathy, lactic acidosis, and stroke-like episodes	Mitochondrial	10	
MT-TS2	Mitochondrial multisystemic disorder	Mitochondrial	2	
MT-TT		Mitochondrial	5	
MT-TV	Hypertrophic cardiomyopathy (HCM), Leigh syndrome, Mitochondrial multisystemic disorder	Mitochondrial	3	
MT-TW	Leigh syndrome, Myopathy, mitochondrial	Mitochondrial	8	
MT-TY	Mitochondrial multisystemic disorder	Mitochondrial	4	
MYH14	Deafness, Peripheral neuropathy, myopathy, hoarseness, and hearing loss	AD	7	44
MYH9	Sebastian syndrome, May-Hegglin anomaly, Epstein syndrome, Fechtner syndrome, Macrothrombocytopenia and progressive sensorineural deafness, Deafness, autosomal dominant 17	AD	25	117

Blueprint Genetics

MYO15A	Deafness	AR	97	235
MYO3A	Deafness	AR	9	22
MYO6	Deafness, Deafness, autosomal dominant, 22	AD/AR	24	68
MYO7A	Deafness, Deafness, autosomal dominant 11, Usher syndrome, type I	AD/AR	239	515
NARS2	Combined oxidative phosphorylation deficiency	AR	12	12
NDP	Exudative vitreoretinopathy, Norrie disease	XL	31	167
NLRP3	Neonatal onset multisystem inflammatory disease (NOMID), Muckle-Wells syndrome, Chronic infantile neurologic cutaneous articular (CINCA) syndrome, Familial cold-induced autoinflammatory syndrome 1	AD	20	136
OSBPL2	Deafness	AD	2	3
<u>OTOA</u> *,#	Deafness	AR	19	28
OTOF	Neuropathy, Deafness	AR	107	163
OTOG	Deafness	AR	18	3
OTOGL	Deafness	AR	26	23
P2RX2	Deafness	AD	2	4
PAX3	Craniofacial-deafness-hand syndrome, Waardenburg syndrome	AD/AR	54	149
PCDH15	Deafness, Usher syndrome, type 1D	AR/Digenic	113	118
PDE1C	Hearing loss	AD	2	2
PDZD7	Deafness, autosomal recessive	AR	11	19
PEX1	Heimler syndrome, Peroxisome biogenesis factor disorder 1A, Peroxisome biogenesis factor disorder 1B	AR	112	134
PEX26	Adrenoleukodystrophy, neonatal, Zellweger syndrome, Peroxisome biogenesis disorder	AR	13	27
PEX6	Heimler syndrome, Peroxisome biogenesis disorder 4A, Peroxisome biogenesis disorder 4B	AR	58	107
PISD		AR		
<u>PNPT1</u> *,#	Deafness, Combined oxidative phosphorylation deficiency, 13	AR	11	13
POLR1C	Treacher Collins syndrome	AR	17	21
POLR1D	Treacher Collins syndrome	AD/AR	9	26
POU3F4	Deafness	XL	25	80
POU4F3	Deafness	AD	9	33
<u>PRPS1</u> *	Phosphoribosylpyrophosphate synthetase I superactivity, Arts syndrome, Charcot-Marie-Tooth disease, X-linked recessive, 5, Deafness, X-linked 1	XL	27	32

Blueprint Genetics

RDX*	Deafness	AR	6	10
RMND1*	Combined oxidative phosphorylation deficiency	AR	17	15
RPS6KA3	Coffin-Lowry syndrome, Mental retardation	XL	65	171
S1PR2	Deafness, autosomal recessive 68	AR	2	3
SALL1*	Townes-Brocks syndrome 1	AD	31	87
SALL4	Acro-renal-ocular syndrome, Duane-radial ray/Okohiro syndrome	AD	21	56
SEMA3E	CHARGE syndrome	AD	1	4
SERPINB6	Deafness	AR	2	3
SIX1	Deafness, Branchiootic syndrome, Branchiootorenal syndrome	AD	11	19
SIX5	Branchiootorenal syndrome	AD	3	10
SLC17A8	Deafness	AD	1	8
SLC19A2	Thiamine-responsive megaloblastic anemia syndrome	AR	14	51
SLC22A4	Hearing loss	AR		2
SLC26A4	Deafness, Pendred syndrome, Enlarged vestibular aqueduct	AR	181	548
SLC26A5	Deafness	AR	2	7
SLC29A3	Histiocytosis-lymphadenopathy plus syndrome, Dysosteosclerosis	AR	17	25
SLC33A1*	Congenital cataracts, hearing loss, and neurodegeneration, Spastic paraplegia 42, autosomal dominant	AD/AR	6	7
SLC52A2	Brown-Vialetto-Van Laere syndrome	AR	27	25
SLC52A3	Fazio-Londe disease, Brown-Vialetto-Van Laere syndrome	AR	30	42
SLITRK6	Deafness and myopia	AR	3	5
SMAD4	Juvenile polyposis/hereditary hemorrhagic telangiectasia syndrome, Polyposis, juvenile intestinal, Myhre dysplasia, Hereditary hemorrhagic telangiectasia	AD	179	143
SMPX	Deafness	XL	8	14
SNAI2	Waardenburg syndrome, Piebaldism	AR	2	4
SOX10	Peripheral demyelinating neuropathy, central dysmyelination, Waardenburg syndrome, and Hirschsprung disease, Kallmann syndrome	AD	56	148
SPATA5	Developmental delay with or without dysmorphic facies and autism, Epilepsy, hearing loss, and mental retardation syndrome (EHLMRS)	AR	27	27
STAG2	Congenital heart defects, dysmorphic facial features, and intellectual developmental disorder	XL	6	14
STRC*,#	Deafness	AR	31	85

Blueprint Genetics

SUCLA2	Mitochondrial DNA depletion syndrome	AR	9	29
SUCLG1	Mitochondrial DNA depletion syndrome	AR	12	28
SYNE4	Deafness	AR	6	2
SYT2	Myasthenic syndrome, congenital 7, presynaptic	AD	3	3
TBC1D24	Deafness, onychodystrophy, osteodystrophy, mental retardation, and seizures (DOORS) syndrome, Deafness, autosomal dominant, 65, Myoclonic epilepsy, infantile, familial, Epileptic encephalopathy, early infantile, 16, Deafness, autosomal recessive 86	AD/AR	43	55
TBL1X	Congenital hypothyroidism, Hearing loss		2	8
TCOF1	Treacher Collins syndrome	AD	50	330
TECTA	Deafness	AD/AR	36	120
TFAP2A	Branchiooculofacial syndrome	AD	23	42
<u>TIMM8A*</u>	Mohr-Tranebjaerg syndrome, Jensen syndrome, Opticoacoustic nerve atrophy with dementia	XL	11	21
TJP2	Cholestasis, progressive familial intrahepatic, Hypercholanemia, familial, Deafness, autosomal dominant 51	AD/AR	25	27
TMC1	Deafness, Deafness, autosomal dominant 36	AD/AR	33	91
TMEM132E	Hearing loss	AR		1
TMIE	Deafness	AR	9	10
TMPRSS3	Deafness	AR	25	82
TNC	Deafness	AD	3	6
TPRN	Deafness	AR	6	12
TRIOBP	Deafness	AR	22	40
TRMU	Liver failure, infantile, Reversible infantile respiratory chain deficiency	AR	20	21
TSHZ1	Aural atresia, congenital	AD	2	4
<u>TSPEAR*</u>	Deafness	AR	2	7
TUBB4B	Leber congenital amaurosis, Hearing loss	AD	2	3
<u>TYR*</u>	Albinism, oculocutaneous	AR	77	441
USH1C	Deafness, Usher syndrome, type IC	AR	45	51
USH1G	Usher syndrome, type 1G	AR	13	32
USH2A	Retinitis pigmentosa 39, Usher syndrome, type 2A	AR	401	1169
VCAN	Wagner disease	AD	11	19
WBP2	Deafness, autosomal recessive 107	AR	3	3



WFS1	Wolfram syndrome, Deafness, Wolfram-like syndrome, autosomal dominant, Deafness, autosomal dominant 6/14/38, Cataract 41	AD/AR	69	362
XYLT2	Spondyloocular syndrome	AR	2	10

*Some regions of the gene are duplicated in the genome. [Read more.](#)

The gene has suboptimal coverage (means <90% of the gene's target nucleotides are covered at >20x with mapping quality score (MQ>20) reads), and/or the gene has exons listed under Test limitations section that are not included in the panel as they are not sufficiently covered with high quality sequence reads.

The sensitivity to detect variants may be limited in genes marked with an asterisk (*) or number sign (#)

Gene refers to the HGNC approved gene symbol; Inheritance refers to inheritance patterns such as autosomal dominant (AD), autosomal recessive (AR), mitochondrial (mi), X-linked (XL), X-linked dominant (XLD) and X-linked recessive (XLR); ClinVar refers to the number of variants in the gene classified as pathogenic or likely pathogenic in this database ([ClinVar](#)); HGMD refers to the number of variants with possible disease association in the gene listed in Human Gene Mutation Database ([HGMD](#)). The list of associated, gene specific phenotypes are generated from [CGD](#) or Mitomap databases.

Non-coding disease causing variants covered by the panel

Gene	Genomic location HG19	HGVS	RefSeq	RS-number
AIFM1	ChrX:129274636	c.697-44T>G	NM_004208.3	
AIFM1	ChrX:129299753	c.-123G>C	NM_004208.3	rs724160014
ANKH	Chr5:14871567	c.-11C>T	NM_054027.4	
BCS1L	Chr2:219524871	c.-147A>G	NM_004328.4	
BCS1L	Chr2:219525123	c.-50+155T>A	NM_004328.4	rs386833855
BTD	Chr3:15683399	c.310-15delT	NM_000060.2	rs587783008
BTD	Chr3:15687154	c.*159G>A	NM_000060.2	rs530872564
CDKN1C	Chr11:2905209	c.*5+20G>T	NM_000076.2	rs760540648
CHD7	Chr8:61734568	c.2836-15C>G	NM_017780.3	
CHD7	Chr8:61757794	c.5051-15T>A	NM_017780.3	
CHD7	Chr8:61763034	c.5405-18C>A	NM_017780.3	rs199981784
CHD7	Chr8:61763035	c.5405-17G>A	NM_017780.3	rs794727423
CHD7	Chr8:61763039	c.5405-13G>A	NM_017780.3	rs1131690787
CLRN1	Chr3:150660197	c.254-649T>G	NM_001195794.1	rs976853535
COL11A1	Chr1:103386637	c.3744+437T>G	NM_080629.2	
COL11A1	Chr1:103488576	c.1027-24A>G	NM_080629.2	
COL11A1	Chr1:103491958	c.781-450T>G	NM_080629.2	rs587782990
COL2A1	Chr12:48379984	c.1527+135G>A	NM_001844.4	



Blueprint Genetics

COL4A3	Chr2:228145145	c.2224-11C>T	NM_000091.4	
COL4A3	Chr2:228168708	c.4028-27A>G	NM_000091.4	
COL4A3	Chr2:228173092	c.4462+457C>G	NM_000091.4	
COL4A3	Chr2:228173596	c.4463-18dupA	NM_000091.4	rs769590145
COL4A4	Chr2:227875240	c.4334-23A>G	NM_000092.4	
COL4A5	ChrX:107813924	c.385-719G>A	NM_033380.2	rs104886396
COL4A5	ChrX:107816792	c.466-12G>A	NM_033380.2	rs104886414
COL4A5	ChrX:107820077	c.609+875G>T	NM_033380.2	
COL4A5	ChrX:107821295	c.646-12_646-11delTT	NM_033380.2	rs104886436
COL4A5	ChrX:107834930	c.1423+57dupC	NM_033380.2	rs104886328
COL4A5	ChrX:107838719	c.1424-20T>A	NM_033380.2	rs281874668
COL4A5	ChrX:107842994	c.1948+894C>G	NM_033380.2	
COL4A5	ChrX:107845097	c.2042-18A>G	NM_033380.2	rs104886341
COL4A5	ChrX:107849932	c.2245-40A>G	NM_033380.2	
COL4A5	ChrX:107849958	c.2245-14T>A	NM_033380.2	
COL4A5	ChrX:107852872	c.2395+2750A>G	NM_033380.2	
COL4A5	ChrX:107908726	c.3374-11C>A	NM_033380.2	rs104886387
COL4A5	ChrX:107933678	c.4529-2300T>G	NM_033380.2	
COL4A5	ChrX:107935633	c.4529-345A>G	NM_033380.2	
COL4A5	ChrX:107938272	c.4821+121T>C	NM_033380.2	rs104886423
COL4A5	ChrX:107938337	c.4822-152dupT	NM_033380.2	
COL4A5	ChrX:107938346	c.4822-151_4822-150insT	NM_033380.2	rs397515494
DCAF17	Chr2:172305176	c.322-14delC	NM_025000.3	rs201494527
DFNA5	Chr7:24746007	c.991-15_991-13delTTC	NM_004403.2	rs727505273
DIAPH3	Chr13:60738072	c.-172G>A	NM_001042517.1	
DIAPH3	Chr13:60738073	c.-173C>T	NM_001042517.1	
EDN3	Chr20:57875743	c.-125G>A	NM_000114.2	
EDN3	Chr20:57875849	c.-19C>A	NM_000114.2	rs375594972
EYA1	Chr8:72156939	c.1051-12T>G	NM_000503.4	
EYA1	Chr8:72211483	c.640-15G>A	NM_000503.4	
EYA4	Chr6:133833847	c.1282-12T>A	NM_004100.4	

EYA4	Chr6:133833997	c.1341-19T>A	NM_004100.4	
GJB2	Chr13:20763744	c.-22-2A>C	NM_004004.5	rs201895089
GJB2	Chr13:20766920	c.-23+2T>A	NM_004004.5	
GJB2	Chr13:20766921	c.-23+1G>A	NM_004004.5	rs80338940
GJB2	Chr13:20766922	c.-23G>T	NM_004004.5	rs786204734
GJB2	Chr13:20767158	c.-259C>T	NM_004004.5	
GJB2	Chr13:20767159	c.-260C>T	NM_004004.5	
GRHL2	Chr8:102505149	c.20+133delA	NM_024915.3	
GRHL2	Chr8:102505272	c.20+257delT	NM_024915.3	
GRHL2	Chr8:102505561	c.20+544G>T	NM_024915.3	
GRXCR1	Chr4:42965170	c.627+19A>T	NM_001080476.2	rs201824235
HGF	Chr7:81384504	c.482+1991_482+2000delGATGATGAAA	NM_000601.4	rs900334407
HGF	Chr7:81384516	c.482+1986_482+1988delTGA	NM_000601.4	
HSD17B4	Chr5:118837725	c.1285-11C>G	NM_001199291.1	rs779466683
KCNJ10	Chr1:160039811	c.-1+1G>T	NM_002241.4	rs796052606
KCNQ1	Chr11:2484803			rs2074238
MYO3A	Chr10:26409593	c.1777-12G>A	NM_017433.4	
MYO6	Chr6:76593963	c.2417-1758T>G	NM_004999.3	
MYO7A	Chr11:76839534	c.-48A>G	NM_000260.3	
MYO7A	Chr11:76893448	c.3109-21G>A	NM_000260.3	
MYO7A	Chr11:76915107	c.5327-14T>G	NM_000260.3	
MYO7A	Chr11:76915110	c.5327-11A>G	NM_000260.3	rs397516316
MYO7A	Chr11:76919448	c.5857-27_5857-26insTTGAG	NM_000260.3	
NDP	ChrX:43818099	c.-207-1G>A	NM_000266.3	
NDP	ChrX:43832545	c.-208+5G>A	NM_000266.3	
NDP	ChrX:43832548	c.-208+2T>G	NM_000266.3	
NDP	ChrX:43832549	c.-208+1G>A	NM_000266.3	
NDP	ChrX:43832685	c.-343A>G	NM_000266.3	rs895911086
NDP	ChrX:43832722	c.-391_-380delCTCTCTCTCCCTinsGTCTCTC	NM_000266.3	
NDP	ChrX:43832724	c.-396_-383delTCCCTCTCTCTCTC	NM_000266.3	rs770996360
PAX3	Chr2:223085913	c.958+28A>T	NM_181459.3	

PCDH15	Chr10:56560684	c.-29+1G>C	NM_001142763.1	
PEX6	Chr6:42933858	c.2301-15C>G	NM_000287.3	rs267608236
PEX6	Chr6:42933952	c.2300+28G>A	NM_000287.3	rs267608237
RPS6KA3	ChrX:20191268	c.1228-279T>G	NM_004586.2	
RPS6KA3	ChrX:20213274	c.326-11A>G	NM_004586.2	
SLC26A4	Chr7:107301201	c.-103T>C	NM_000441.1	rs60284988
SLC26A4	Chr7:107301244	c.-60A>G	NM_000441.1	rs545973091
SLC26A4	Chr7:107301301	c.-4+1G>C	NM_000441.1	
SLC26A4	Chr7:107301305	c.-4+5G>A	NM_000441.1	rs727503425
SLC26A4	Chr7:107323842	c.918+45_918+47delCAA	NM_000441.1	
SLC26A4	Chr7:107330533	c.1150-35_1150-28delTTTGTAGG	NM_000441.1	
SLC26A4	Chr7:107334836	c.1264-12T>A	NM_000441.1	
SLC26A4	Chr7:107336364	c.1438-7dupT	NM_000441.1	rs754734032
SLC26A4	Chr7:107341513	c.1708-27_1708-11delTAAGTAACTTGACATTT	NM_000441.1	
SLC26A4	Chr7:107350439	c.2090-52_2090-49delCAAA	NM_000441.1	
SLC29A3	Chr10:73122778	c.*413G>A	NM_018344.5	
SLC52A2	Chr8:145582843	c.-110-1G>A	NM_024531.4	
SNAI2	Chr8:49833972	c.-149_-148delCGinsTA	NM_003068.4	
SOX10	Chr22:38379877	c.-84-2A>T	NM_006941.3	
SOX10	Chr22:38412215	c.-31954C>T	NM_006941.3	rs606231342
SOX10	Chr22:38412781	c.-32520C>G	NM_006941.3	rs533778281
TIMM8A	ChrX:100601671	c.133-23A>C	NM_004085.3	rs869320666
TMC1	Chr9:75315577	c.362+18A>G	NM_138691.2	
TYR	Chr11:88960973	c.1037-18T>G	NM_000372.4	
USH2A	Chr1:215821092	c.14583-20C>G	NM_206933.2	
USH2A	Chr1:215967783	c.9959-4159A>G	NM_206933.2	
USH2A	Chr1:216039721	c.8845+628C>T	NM_206933.2	
USH2A	Chr1:216064540	c.7595-2144A>G	NM_206933.2	rs786200928
USH2A	Chr1:216247476	c.5573-834A>G	NM_206933.2	
USH2A	Chr1:216592035	c.486-14G>A	NM_206933.2	rs374536346
USH2A	Chr1:216596610	c.-259G>T	NM_206933.2	

Test Strengths

Our panel assay enables the detection of common deletions in GJB6 such as (~309 kb del (GJB6-D13S1830) and ~232 kb del (GJB6-D13S1854)).

The strengths of this test include:

- CAP accredited laboratory
- CLIA-certified personnel performing clinical testing in a CLIA-certified laboratory
- Powerful sequencing technologies, advanced target enrichment methods and precision bioinformatics pipelines ensure superior analytical performance
- Careful construction of clinically effective and scientifically justified gene panels
- Some of the panels include the whole mitochondrial genome (please see the Panel Content section)
- Our Nucleus online portal providing transparent and easy access to quality and performance data at the patient level
- Our publicly available analytic validation demonstrating complete details of test performance
- ~2,000 non-coding disease causing variants in our clinical grade NGS assay for panels (please see 'Non-coding disease causing variants covered by this panel' in the Panel Content section)
- Our rigorous variant classification scheme
- Our systematic clinical interpretation workflow using proprietary software enabling accurate and traceable processing of NGS data
- Our comprehensive clinical statements

Test Limitations

Variants in the *KCNE1* gene should not be used for risk assessment at the moment. Specifically, *KCNE1* c.253G>A, p.(Asp85Asn) variant has been considered to be a mild risk factor for acquired long QT syndrome. However, in the newest version of the reference genome GRCh38, a gene *KCNE1B*, nearly identical to *KCNE1* has appeared. By using standard NGS technologies, as well as Sanger sequencing, it is not possible to get reliable region-specific sequences for these genes. It is likely that reads that have been earlier mapped to *KCNE1* actually belong to *KCNE1B*. Moreover, it is currently unclear whether *KCNE1B* produces a protein product, and if a protein is produced, whether the gene is expressed in heart. More independent data characterizing *KCNE1B* and its function are needed. Currently, all *KCNE1* sequence data and the literature related to *KCNE1* variants should be interpreted with caution.

The following exons are not included in the panel as they are not sufficiently covered with high quality sequence reads: *OTOA* (NM_144672:22-27), *PDZD7* (NM_024895:10), *POLR1C* (NM_001318876:9), *STRC* (NM_153700:1-18). Note that we are able to detect variant in exons 19-29 of *STRC*, but our abilities are limited due to the high degree of homology that is shared between these exons and other regions of the genome. Whole gene deletions of *STRC* can and have been detected.

Genes with suboptimal coverage in our assay are marked with number sign (#) and genes with partial, or whole gene, segmental duplications in the human genome are marked with an asterisk (*) if they overlap with the UCSC pseudogene regions. Gene is considered to have suboptimal coverage when >90% of the gene's target nucleotides are not covered at >20x with mapping quality score (MQ>20) reads. The technology may have limited sensitivity to detect variants in genes marked with these symbols (please see the Panel content table above).

This test does not detect the following:

- Complex inversions
- Gene conversions
- Balanced translocations
- Some of the panels include the whole mitochondrial genome but not all (please see the Panel Content section)
- Repeat expansion disorders unless specifically mentioned

Blueprint Genetics

- Non-coding variants deeper than ± 20 base pairs from exon-intron boundary unless otherwise indicated (please see above Panel Content / non-coding variants covered by the panel).

This test may not reliably detect the following:

- Low level mosaicism in nuclear genes (variant with a minor allele fraction of 14.6% is detected with 90% probability)
- Stretches of mononucleotide repeats
- Low level heteroplasmy in mtDNA (>90% are detected at 5% level)
- Indels larger than 50bp
- Single exon deletions or duplications
- Variants within pseudogene regions/duplicated segments
- Some disease causing variants present in mtDNA are not detectable from blood, thus post-mitotic tissue such as skeletal muscle may be required for establishing molecular diagnosis.

The sensitivity of this test may be reduced if DNA is extracted by a laboratory other than Blueprint Genetics.

For additional information, please refer to the Test performance section and see our Analytic Validation.

Test performance

The Blueprint Genetics comprehensive hearing loss and deafness panel covers classical genes associated with Waardenburg syndrome, Alport syndrome, sensorineural hearing loss, unilateral and bilateral, non-syndromic genetic deafness, Pendred syndrome, Usher syndrome, Stickler syndrome, Jervell and Lange-Nielsen syndrome, Mohr-Tranebjaerg syndrome, Norrie disease, Treacher Collins syndrome, CHARGE syndrome and Branchio-oto-renal (BOR) syndrome. The genes on the panel have been carefully selected based on scientific literature, mutation databases and our experience.

Our panels are sliced from our high-quality whole exome sequencing data. Please see our sequencing and detection performance table for different types of alterations at the whole exome level (Table).

Assays have been validated for different starting materials including EDTA-blood, isolated DNA (no FFPE), saliva and dry blood spots (filter card) and all provide high-quality results. The diagnostic yield varies substantially depending on the assay used, referring healthcare professional, hospital and country. Blueprint Genetics' Plus Analysis (Seq+Del/Dup) maximizes the chance to find a molecular genetic diagnosis for your patient although Sequence Analysis or Del/Dup Analysis may be a cost-effective first line test if your patient's phenotype is suggestive of a specific mutation type.

The genes on the panel have been carefully selected based on scientific literature, mutation databases and our experience.

Our panels are sectioned from our high-quality, clinical grade NGS assay. Please see our sequencing and detection performance table for details regarding our ability to detect different types of alterations (Table).

	Sensitivity % (TP/(TP+FN))	Specificity %
Single nucleotide variants	99.89% (99,153/99,266)	>99.9999%
Insertions, deletions and indels by sequence analysis		
1-10 bps	99.2% (7,745/7,806)	>99.9999%
11-50 bps	99.13% (2,524/2,546)	>99.9999%
Copy number variants (exon level dels/dups)		
1 exon level deletion (heterozygous)	100% (20/20)	NA
1 exon level deletion (homozygous)	100% (5/5)	NA
1 exon level deletion (het or homo)	100% (25/25)	NA

Blueprint Genetics



2-7 exon level deletion (het or homo)	100% (44/44)	NA
1-9 exon level duplication (het or homo)	75% (6/8)	NA
Simulated CNV detection		
5 exons level deletion/duplication	98.7%	100.00%
Microdeletion/-duplication sdrs (large CNVs, n=37)		
Size range (0.1-47 Mb)	100% (25/25)	

The performance presented above reached by Blueprint Genetics high-quality, clinical grade NGS sequencing assay with the following coverage metrics

Mean sequencing depth	143X
Nucleotides with >20x sequencing coverage (%)	99.86%

Performance of Blueprint Genetics Mitochondrial Sequencing Assay.

	Sensitivity %	Specificity %
ANALYTIC VALIDATION (NA samples; n=4)		
Single nucleotide variants		
Heteroplasmic (45-100%)	100.0% (50/50)	100.0%
Heteroplasmic (35-45%)	100.0% (87/87)	100.0%
Heteroplasmic (25-35%)	100.0% (73/73)	100.0%
Heteroplasmic (15-25%)	100.0% (77/77)	100.0%
Heteroplasmic (10-15%)	100.0% (74/74)	100.0%
Heteroplasmic (5-10%)	100.0% (3/3)	100.0%
Heteroplasmic (<5%)	50.0% (2/4)	100.0%
CLINICAL VALIDATION (n=76 samples)		
All types		
Single nucleotide variants n=2026 SNVs		
Heteroplasmic (45-100%)	100.0% (1940/1940)	100.0%
Heteroplasmic (35-45%)	100.0% (4/4)	100.0%
Heteroplasmic (25-35%)	100.0% (3/3)	100.0%
Heteroplasmic (15-25%)	100.0% (3/3)	100.0%
Heteroplasmic (10-15%)	100.0% (9/9)	100.0%



Blueprint Genetics



Heteroplasmic (5-10%)	92.3% (12/13)	99.98%
Heteroplasmic (<5%)	88.9% (48/54)	99.93%
Insertions and deletions by sequence analysis n=40 indels		
Heteroplasmic (45-100%) 1-10bp	100.0% (32/32)	100.0%
Heteroplasmic (5-45%) 1-10bp	100.0% (3/3)	100.0%
Heteroplasmic (<5%) 1-10bp	100.0% (5/5)	99,997%
SIMULATION DATA /(mitomap mutations)		
Insertions, and deletions 1-24 bps by sequence analysis; n=17		
Homoplasmic (100%) 1-24bp	100.0% (17/17)	99.98%
Heteroplasmic (50%)	100.0% (17/17)	99.99%
Heteroplasmic (25%)	100.0% (17/17)	100.0%
Heteroplasmic (20%)	100.0% (17/17)	100.0%
Heteroplasmic (15%)	100.0% (17/17)	100.0%
Heteroplasmic (10%)	94.1% (16/17)	100.0%
Heteroplasmic (5%)	94.1% (16/17)	100.0%
Copy number variants (separate artificial mutations; n=1500)		
Homoplasmic (100%) 500 bp, 1kb, 5 kb	100.0%	100.0%
Heteroplasmic (50%) 500 bp, 1kb, 5 kb	100.0%	100.0%
Heteroplasmic (30%) 500 bp, 1kb, 5 kb	100.0%	100.0%
Heteroplasmic (20%) 500 bp, 1kb, 5 kb	99.7%	100.0%
Heteroplasmic (10%) 500 bp, 1kb, 5 kb	99.0%	100.0%
The performance presented above reached by following coverage metrics at assay level (n=66)		
	Mean of medians	Median of medians
Mean sequencing depth MQ0 (clinical)	18224X	17366X
Nucleotides with >1000x MQ0 sequencing coverage (%) (clinical)	100%	
rho zero cell line (=no mtDNA), mean sequencing depth	12X	

Bioinformatics

The target region for each gene includes coding exons and ± 20 base pairs from the exon-intron boundary. In addition, the panel includes non-coding variants if listed above (Non-coding variants covered by the panel). Some regions of the gene(s) may be removed from the panel if specifically mentioned in the 'Test limitations' section above. The sequencing data generated in our laboratory is analyzed with our proprietary data analysis and annotation pipeline, integrating state-of-the-art



Blueprint Genetics

algorithms and industry-standard software solutions. Incorporation of rigorous quality control steps throughout the workflow of the pipeline ensures the consistency, validity and accuracy of results. Our pipeline is streamlined to maximize sensitivity without sacrificing specificity. We have incorporated a number of reference population databases and mutation databases such as, but not limited to, [1000 Genomes Project](#), [gnomAD](#), [ClinVar](#) and [HGMD](#) into our [clinical interpretation software to make the process effective and efficient](#). For missense variants, *in silico* variant prediction tools such as SIFT, PolyPhen, MutationTaster are used to assist with variant classification. Through our [online ordering and statement reporting system, Nucleus](#), the customer has an access to details of the analysis, including [patient specific sequencing metrics, a gene level coverage plot and a list of regions with inadequate coverage if present](#). This reflects our mission to [build fully transparent diagnostics where customers have easy access to crucial details of the analysis process](#).

Clinical interpretation

We provide customers with the most comprehensive clinical report available on the market. Clinical interpretation requires a fundamental understanding of clinical genetics and genetic principles. At Blueprint Genetics, our PhD molecular geneticists, medical geneticists and clinical consultants prepare the clinical statement together by evaluating the identified variants in the context of the phenotypic information provided in the requisition form. Our goal is to provide clinically meaningful statements that are understandable for all medical professionals regardless of whether they have formal training in genetics.

Variant classification is the corner stone of clinical interpretation and resulting patient management decisions. Our classifications follow the [ACMG guideline 2015](#).

The final step in the analysis of sequence variants is confirmation of variants classified as pathogenic or likely pathogenic using bi-directional Sanger sequencing. Variant(s) fulfilling the following criteria are not Sanger confirmed: the variant quality score is above the internal threshold for a true positive call, and visual check-up of the variant at IGV is in-line with the variant call. Reported variants of uncertain significance are confirmed with bi-directional Sanger sequencing only if the quality score is below our internally defined quality score for true positive call. Reported copy number variations with a size <10 exons are confirmed by orthogonal methods such as qPCR if the specific CNV has been seen less than three times at Blueprint Genetics.

Our clinical statement includes tables for sequencing and copy number variants that include basic variant information (genomic coordinates, HGVS nomenclature, zygosity, allele frequencies, *in silico* predictions, OMIM phenotypes and classification of the variant). In addition, the statement includes detailed descriptions of the variant, gene and phenotype(s) including the role of the specific gene in human disease, the mutation profile, information about the gene's variation in population cohorts and detailed information about related phenotypes. We also provide links to the references used, congress abstracts and mutation variant databases used to help our customers further evaluate the reported findings if desired. The conclusion summarizes all of the existing information and provides our rationale for the classification of the variant.

Identification of pathogenic or likely pathogenic variants in dominant disorders or their combinations in different alleles in recessive disorders are considered molecular confirmation of the clinical diagnosis. In these cases, family member testing can be used for risk stratification within the family. In the case of variants of uncertain significance (VUS), we do not recommend family member risk stratification based on the VUS result. Furthermore, in the case of VUS, we do not recommend the use of genetic information in patient management or genetic counseling.

Our interpretation team analyzes millions of variants from thousands of individuals with rare diseases. Thus, our database, and our understanding of variants and related phenotypes, is growing by leaps and bounds. Our laboratory is therefore well positioned to re-classify previously reported variants as new information becomes available. If a variant previously reported by Blueprint Genetics is re-classified, our laboratory will issue a follow-up statement to the original ordering health care provider at no additional cost.

#}

ICD codes

Commonly used ICD-10 codes when ordering the Comprehensive Hearing Loss and Deafness Panel

ICD-10	Disease
--------	---------

Blueprint Genetics



E70.30	Waardenburg syndrome
Q87.89	Alport syndrome
F84.2	Rett syndrome
H90.5	Sensorineural hearing loss, unilateral and bilateral
E07.1	Pendred syndrome
H35.50	Usher syndrome, type IV
Q89.8	Stickler syndrome
G31.89	Mohr-Tranebjaerg syndrome
H35.50	Norrie disease
Q75.4	Treacher Collins syndrome
Q89.8	CHARGE syndrome
H49.40	Progressive external ophthalmoplegia
Q87.89	Branchio-oto-renal (BOR) syndrome
G11.9	Hereditary ataxia
C94.2	Acute Megakaryoblastic Leukemia
K59.8	Chronic Intestinal Pseudoobstruction
T36.5	Adverse effect of aminoglycosides
G93.41	Metabolic Encephalopathy
H49.81	Kearns Sayre Syndrome
E88.42	MERFF Syndrome
H47.013	Nonarteritic Anterior Ischemic Optic Neuropathy
G60.2	Neuropathy in association with hereditary ataxia
G30	Alzheimer's Disease
G25.5	Chorea
G40	Epilepsy and recurrent seizures
I42	Cardiomyopathy
N26.9	Focal Segmental Glomerulosclerosis
G31.82	Leigh's Disease
H47.2	Leber's hereditary optic neuropathy
G71.3	Mitochondrial Myopathy
I42.1	Hypertrophic Cardiomyopathy



E11.9	Non-Insulin Dependent Diabetes Mellitus
Z86.74	Personal history of sudden cardiac arrest
H90.3	Sensorineural Hearing Loss

Accepted sample types

- EDTA blood, min. 1 ml
- Purified DNA, min. 3µg*
- Saliva (Oragene DNA OG-500 kit)

Label the sample tube with your patient's name, date of birth and the date of sample collection.

Note that we do not accept DNA samples isolated from formalin-fixed paraffin-embedded (FFPE) tissue.

Resources

- [Alford RL et al. American College of Medical Genetics and Genomics guideline for the clinical evaluation and etiologic diagnosis of hearing loss. Genet Med. 2014 Apr;16\(4\):347-55.](#)
- [Alport Syndrome Foundation](#)
- [American Hearing Research Foundation](#)
- [American Society for Deaf Children](#)
- [British Deaf Association](#)
- [CHARGE Syndrome Foundation](#)
- [GeneReviews - BOR Syndrome](#)
- [GeneReviews - Deafness and Hereditary Hearing Loss](#)
- [GeneReviews - Pendred Syndrome](#)
- [GeneReviews - Stickler Syndrome](#)
- [GeneReviews - Treacher Collins syndrome](#)
- [GeneReviews - Usher Syndrome, type I](#)
- [GeneReviews - Usher Syndrome, type II](#)
- [GeneReviews - Waardenburg Syndrome](#)
- [Hearing Health Foundation](#)
- [Hearing Link](#)
- [Hereditary Hearing loss](#)
- [NORD - Alport Syndrome](#)
- [NORD - BOR Syndrome](#)
- [NORD - CHARGE Syndrome](#)
- [NORD - Jervell and Lange-Nielsen Syndrome](#)
- [NORD - Norrie Disease](#)
- [NORD - Stickler Syndrome](#)
- [NORD - Treacher Collins syndrome](#)
- [NORD - Usher Syndrome](#)
- [NORD - Waardenburg Syndrome](#)
- [National Association of the Deaf US](#)
- [Norrie Disease Association](#)
- [Sloan-Heggen CM et al. Navigating genetic diagnostics in patients with hearing loss. Curr Opin Pediatr. 2016 Dec;28\(6\):705-712.](#)
- [Stickler Syndrome Support Group UK](#)
- [Stickler Syndrome Support Group US](#)
- [The Cain Foundation for BOR syndrome](#)
- [Usher Syndrome Coalition](#)